




Sizing Guide

TADAWI
LIFE TECH

Your Life...Our Commitment

حياتكم.. التزامنا

 Muscat, Oman

 +968-91235160, 79012020


 info@tadawilifetech.com

 Tadawilifetech.com



The sizing of all Ortino back braces (**LO**, **LSO**, **TLSO**, **TLSO Plus**) is determined based on the waist circumference.


Sizes	Back height (cm)
M	41 – 46
L	47 – 53
XL	54 – 60
2XL	66 – 61
3XL	66 – 61





For the **Spinomed back brace**, in addition to the waist circumference, the height of the spine must also be measured. Therefore, the sizing of the Spinomed is based on two measurements.


Sizes	Back height (cm)
S	41 – 46
M	47 – 53
L	54 – 60
XL	66 – 61





- The sizing of the **Wrist Brace** is determined based on the wrist circumference.

Sizes	Wrist circumference (cm)
S	12 – 15.5
M	16 – 17.5
L	18 – 23





- The sizing of the **Tennis Elbow Brace** is based on the forearm circumference.


Sizes	Forearm circumference (cm)
Size 1	20 – 28
Size 2	29 – 32





- The sizing of the **Rehabilitation knee brace** is based on the thigh circumference measured 15 cm above the knee.


Sizes	Thigh circumference (cm)
M	* < 44
L	45 – 57
XL	* > 58





- The sizing of the **Rehabilitation knee brace** is based on the thigh circumference measured 15 cm above the knee.


Sizes	Thigh circumference (cm)
S	42 – 47
M	48 – 52
L	53 – 57
XL	62 – 58
2XL	63 – 69





- The sizing of the **Patella Brace** is based on the knee circumference, measured just below the knee.

Sizes	The circumference measurement below the knee (cm)
S	27 – 32
M	33 – 38
L	39 – 43





- The sizing of the **Ankle Brace** is based on the ankle circumference.

Sizes	Ankle Circumference (cm)
M	19 -22
L	23- 28
XL	* > 29





- The sizing of the **Foot Drop Brace** is also based on the ankle circumference.

Sizes	Ankle Circumference (cm)
Size 1	19 – 24
Size 2	25 – 32





● The following products are **Free Size**.

- Post-Op Knee Brace
- Unloader Knee Brace
- Immobilizer Knee Brace
- Posture Corrector
- Post-Op Elbow Brace (Left & Right)